




Phone: 954-772-5052

Email: LSOLOGUREN@AHCEQUIP.COM



# Compliance Report


Initial compliance period	05/04/2023 - 06/02/2023
Compliance met	Yes
Compliance percentage	93%
Payor	Medicare

Usage	05/04/2023 - 06/02/2023
Usage days	29/30 days (97%)
>= 4 hours	28 days (93%)
< 4 hours	1 days (3%)
Usage hours	210 hours 21 minutes
Average usage (total days)	7 hours 1 minutes
Average usage (days used)	7 hours 15 minutes
Median usage (days used)	7 hours 21 minutes
Total used hours (value since last reset - 06/02/2023)	213 hours

AirSense 11 AutoSet	
Serial number	23223608894
Mode	AutoSet
Min Pressure	5 cmH2O
Max Pressure	20 cmH2O
EPR	Ramp Only
EPR level	3
Response	Standard

Therapy						
Pressure - cmH2O	Median:	8.9	95th percentile: 1	.3.3	Maximum:	15.5
Leaks - L/min	Median:	13.8	95th percentile: 4	12.6	Maximum:	75.4
Events per hour	AI:	1.7	HI:	0.4	AHI:	2.1
Apnea Index	Central:	0.0	Obstructive:	1.5	Unknown:	0.1
RERA Index						0.2
Cheyne-Stokes respiration (average duration per night)			0 minutes	s (0%)		

## Usage - hours





Phone: 954-772-5052


Email: LSOLOGUREN@AHCEQUIP.COM




## Therapy Report

AirSense 11 AutoSet SN: 23223608894


### Usage (hours)




#### Leak (L/min)



#### Pressure (cmH2O)



### AHI (events/hour)





Phone: 954-772-5052

Email: LSOLOGUREN@AHCEQUIP.COM




## Therapy Report

AirSense 11 AutoSet SN: 23223608894


#### Tidal Volume(ml)

Maximum (avg) 95th % (avg) Median (avg)




### Respiratory Rat@breaths/min)

Maximum (avg) 95th % (avg) Median (avg)



## Minute Ventilation(L/min)

Maximum (avg) 95th % (avg) Median (avg)








Phone: 954-772-5052

Email: LSOLOGUREN@AHCEQUIP.COM



SN: 23223608894

## Therapy Report

AirSense 11 AutoSet

### Care Check-In

No available data because the user declined to answer the therapy questions, this feature was turned off on the user's machine or the user did not consent to share their data with myAir.

#### About Care Check-In

Short patient surveys are conducted on the AirSense 11 and myAir app to better understand how the patient is doing with therapy. If the patient consents, we can share Care Check-In information through AirView for clinical users.

#### Questions and answers

At the beginning of therapy, the patient is asked a baseline question. A few additional questions are asked, about once per week, during the early stages of therapy. For example:

How sleepy did you feel this week? Extremely, Very, Moderately, Slightly, Not at all.

How is your therapy? Great, OK, Challenging.

What bothers you most? Getting used to therapy, Fitting the mask, Using the machine.